

MICROLUX/DL SPECIFICATIONS

INTENDED USE: To be used in combination with traditional oral examination by incandescent light by a dentist or health care provider to increase identification, evaluation, and monitoring of oral mucosal abnormalities.

LENGTH: 5" X 5/8" diameter (13cm x 1.5cm diameter)

WEIGHT: 2oz (57 grams)

BATTERY TYPE: N cells (3) Carbon Zinc or Alkaline

BATTERY LIFE: Approx. 40 hours


LIGHT SOURCE: T-1^{3/4} White LED

LIGHT GUIDE: Diffused Glass Fiber Optic

MOUTH RINSE: Acetic Acid (30ml), 1% dilute solution



Acetic Acid
(Mouth rinse)



Light Guide

Nose Cone

Battery Cap &
On/Off Switch