

MICROLUX TRANSILLUMINATOR SPECIFICATIONS

INTENDED USE: To detect cracks, caries, sub gingival calculus and root canal orifice

Length: 5' x 5/8" diameter (13cm x 1.5cm diameter)


Weight: 2.0 oz (57 grams)

Output: 90mw/cm² @ 5600K

Battery: N cells (3) Carbon Zinc or Alkaline. Approx. 40 hrs operation

Light guide: Fiber Optic with turbo action, exit diameter .120 inches

Light Source: T-1^{3/4} White LED



Light guide

Nose Cone

Battery Cap & On/Off Switch