

**MICROLUX LIGHTED MIRROR  
INSTRUCTION LEAFLET**

The Microlux Lighted Mirror is a fiber optic attachment for the Microlux Transilluminator. The Lighted Mirror provides additional light for hard to reach areas. It has a front surface mirror design for clearer vision.


Insert the mirror attachment into the Microlux using firm pressure.

**Maintenance & Cleaning:**

The Lighted Mirror can be autoclaved or cleaned with disinfectant solution. A Chemclave is not recommended.



AdDent, Inc • 43 Miry Brook Road • Danbury, CT 06810 • Tel: 203-778-0200 • Fax: 203-792-2275



**EMERGO EUROPE**  
Molenstraat 15  
2513 BH, The Hague  
The Netherlands

