

## TRIMAX COMPOSITE INSTRUMENT SPECIFICATIONS

**INTENDED USE:** For posterior composite placement

### **HANDLE**

Length: 5.5" x 3/8" (140 x 10 mm) – High Temperature Polysulfone Plastic

Weight: 0.5 oz (14 grams)

### **Kit**

Length: 7" x 4" x 1" (178 x 102 x 25 mm)

Weight: 5.6 oz (158 grams)

### **Tips**

3 Sizes – Polycarbonate Plastic

