

CoMax composite dispenser by AdDent

I have always been an advocate for investing in quality for products you use every day. That might be the perfect pair of comfortable, stylish dress shoes, the perfect navy sports coat, or even a great television set that you settle down in front of in the evenings. I have no problem investing in things I use every day. The same goes for the dental practice. There are things we use daily, and there are things we use maybe just a few times each year. Find bargains on the products you use sporadically. For those that you use daily, get stuff you like.

Recently I found myself annoyed with many of the composite guns in my practice. Most of them were freebies from dental meetings. All of them were made of cheap-feeling, flimsy plastic and did not securely hold a compule of composite. After being annoyed and frustrated somewhere around 1,000 times, I decided to invest in a quality composite dispenser. The CoMax from AdDent was the one I chose.

The CoMax composite dispenser is the polar opposite of the crappy plastic giveaway dispensers. The first thing you will notice is how solid the CoMax feels in your hands. It has a comfortable, reassuring weight to it. This is not to say that it feels heavy—but more like it is well constructed and durable.



One of my main annoyances with all the other composite guns in my office was how easily a compule of composite would fall out or lodge itself in an unusable position. The plunger on most plastic composite guns doesn't actually hold the composite in place. On the CoMax, the compule is held firmly in place with no way to wiggle or fall out. With a quick pull on the neck, the spent compule can be removed and replaced with a fresh one. The mechanism for

this is well designed, easy to use, and secure. It can even be done with one hand.

The trigger on the CoMax is designed to allow for easier dispensing and is very ergonomic. It also allows you to fully express the entire contents of a compule. Many of the lower-quality dispensers tend to bottom out and leave unused material in the compule. Composite is expensive! I don't love the idea of throwing a lot of it away.

This dispenser is fully

autoclavable, and quite frankly, it can take a beating. Instead of having 10 plastic dispensers, we have four CoMax dispensers. They have held up immaculately in spite of heavy use, and they show zero signs of wear and tear.

I didn't know how much I needed a quality composite dispensing gun. It was not until I used a CoMax that I knew what a difference it could make. Give a CoMax a try and you won't go back. Line-drive single to left field for AdDent!