

# CoMax

COMPOSITE DISPENSER



43 Miry Brook Road • Danbury, CT 06810 • (203) 778-0200

**AdDent**  
INCORPORATED

Advancing Dental  
Excellence

ISO 9001 CERTIFIED & 13485 COMPLIANT | ALL PRODUCTS ARE FDA CLEARED

Made in the USA | [www.addent.com](http://www.addent.com)

## COMAX COMPOSITE DISPENSER


The CoMax has a special patented ergonomic mechanism for comfortable handling. The handle minimizes the dispensing force needed to express heavily filled composites. The top flange is pulled back to insert a compule, and easily pulled back again to release it. The CoMax accepts most mfg.'s composite compules. The unit is autoclaveable, and carries a 2 year warranty. The CoMax Composite Dispenser received a 5 Star Award from Reality Publishing.

### FEATURES

- Accepts most composite compules
- Ergonomic handle and trigger
- High mechanical advantage mechanism
- Easy one hand unloading
- Hi tech plastic and stainless steel materials
- Completely autoclaveable
- Expresses all Composite
- Warmed composite provide quicker cures
- Less Microleakage
- Stronger restoration
- Less shrinkage stress

### BENEFITS

- Universal operation
- Easy to hold and use
- Lower dispensing force required
- Time saving operation
- Robust long life operation
- Improved asepsis
- Saves money



Loading the compule into the CoMax dispenser



Holding the CoMax to express composite



Removing the used compule from the CoMax dispenser

Follow us on social media: [f](#) [in](#) [t](#) [v](#) [info@addent.com](mailto:info@addent.com)



**AdDent**  
INCORPORATED

Advancing Dental  
Excellence

ISO 9001 CERTIFIED & 13485 COMPLIANT | ALL PRODUCTS ARE FDA CLEARED



Made in the USA | [www.addent.com](http://www.addent.com)